IPS Empress CAD

[en] Instructions for Use

Leucite glass-ceramic block

[de] Gebrauchsinformation

Leuzit-Glaskeramik-Block

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Patients with permanent teeth

Intended users

- Dentists (chairside fabrication of restorations; clinical workflow)
- Dental laboratory technicians (fabrication of restorations in the dental laboratory)

Special training

No special training required.

Use

For dental use only.

Description

IPS Empress® CAD is a tried-and-tested leucite-reinforced glass-ceramic block, consisting of a glass and a crystal phase, for the fabrication of fixed, full-contour single-tooth restorations in anterior and posterior teeth. IPS Empress CAD can be milled in an authorized CAD/CAM^[1] machine. After wet-milling of the block, the restoration is polished with e.g. OptraGloss® or individually characterized and glazed. For individual characterization and glazing, use IPS e.max® Crystall./ materials "chairside" or IPS Ivocolor materials "labside". The mean biaxial flexural strength from 10 years of quality measurements is 185 MPa (according to ISO 6872:2015).

Technical specifications

Property	Specifications	Typical average value
CTE (25–500°C) [10 ⁻⁶ K ⁻¹]	175 ± 0.5	–
Flexural strength (biaxial) [MPa]	≥ 140	185
Chemical solubility [µg/cm ²]	< 100	–
Type / Class	Type II / Class 2	–

According to ISO 6872:2024

Indications

- Missing tooth structure in anterior and posterior teeth

Types of restorations:

- Veneers
- Inlays
- Onlays (e.g. partial crowns)
- Crowns

Contraindications

The use of the product is contraindicated if the patient is known to be allergic to any of its ingredients.

- Patients with substantially reduced residual dentition

Limitations of use

- Veneered anterior and posterior crowns
- Conventional and self-adhesive cementation
- Temporary cementation
- Very deep sub-gingival preparations
- Untreated bruxism (a splint is indicated after incorporation)
-  Do not re-use.

In the following situations, a successful procedure cannot be ensured:

- Falling short of the required minimum thicknesses
- Machining of the blocks in a non-compatible CAD/CAM system^[1]
- IPS Empress CAD restorations must not be stained and glazed using other materials than the approved ones.
- IPS e.max CAD Crystall./Shades, Stains, Glaze and IPS Ivocolor Shades, Essence, Glaze must not be mixed with one another or applied one after the other.
- When using IPS e.max CAD Crystall./Add-On

IPS e.max CAD Crystall./Glaze Spray:

- Intraoral application
- Restorations whose inner surfaces cannot be reliably and precisely covered (e.g. inlays)

Side effects

There are no known side effects to date.

Interactions

There are no known interactions to date.

Clinical benefit

- Reconstruction of chewing function
- Restoration of esthetics

Composition

Dental ceramics (leucite glass-ceramics), inorganic pigments

2 Application

Block concept

All blocks are made of the same material and feature the same physical properties. Due to esthetic reasons, the individual blocks (translucency levels) are recommended for the following types of restorations:

Degree of translucency	Block sizes/ shades	Processing technique		Types of restorations			
		Polishing technique 	Staining technique 	Veneer	Inlay, onlay	Partial crown	Anterior and Posterior crown
Multi	I12, C14, C14L; shades BL1, BL3, A1, A2, A3, A3.5, B1 ^(a)	✓	✓	✓		✓	✓
HT High Translucency	I10, I12; Shades A1, A2, A3, A3.5, B1, B2, C2, D3 ^(a)	✓	✓	✓	✓		
LT Low Translucency	I12, C14; shades BL1, BL2, BL3, BL4, A1, A2, A3, A3.5, B1, B2, B3, C2, D3 ^(a)	✓	✓	✓		✓	✓

^(a) The availability of shades may vary depending on the degree of translucency/block size or CAD/CAM machine.

Compatible luting materials

Adhesive cementation (e.g. Variolink® Esthetic, Multilink® Automix)

For further information, please refer to the Flow Chart "Pretreatment and cementation".


Shade selection

Clean the teeth prior to shade determination. Determine the shade based on the moist and unprepared tooth or the neighbouring teeth. Determine the shade of the prepared tooth based on the largest discoloured area of the preparation.

Preparation

Prepare the tooth based on the guidelines for all-ceramic restorations and make sure to observe the minimum wall thicknesses:

- No angles or sharp edges
- Shoulder preparation with rounded inner edges and/or chamfer preparation
- The indicated dimensions reflect the minimum thickness for IPS Empress CAD restorations.
- The thickness of the incisal edge of the preparation, particularly in anterior crowns, must be at least 1.0 mm in order to ensure optimum processing in the CAD/CAM machine.

Dimensions in mm

Finishing

Observe the grinding instrument recommendations^[2] and the minimum layer thicknesses when finishing and contouring the ceramic structures. Overheating of the ceramic must be avoided. Smooth out the attachment point of the block, paying particular attention to the proximal contacts. If required, carry out individual shape adjustments. Finish the functional areas (occlusal contact surfaces) of the restoration with a fine diamond to smooth out the surface structure created by the CAD/CAM process. Make sure to thoroughly clean the restoration before further processing and to remove any residue of the milling additive of the CAD/CAM milling unit. Residue of the milling additive remaining on the surface may result in bonding problems and discolouration. If desired, the restoration can be carefully tried in intraorally in order to adjust the occlusion/articulation. Afterwards, the restoration should be cleaned again. The restoration must not be blasted with Al_2O_3 or glass polishing beads. When the restoration is glazed, it may be advisable to create a more pronounced surface texture as glazing will smoothen it.

Completion

– Polishing technique

OptraGloss and/or dental lab polishers can be used for extraoral polishing. Please observe the grinding instrument recommendations^[2]. Avoid overheating of the restoration. Polish proximal areas and large surfaces extraorally prior to cementation. After having cemented the restoration, adjust occlusion/articulation using suitable grinding instruments^[2], e.g. OptraGloss. Subsequently, polish the restoration according to the grinding instrument recommendations^[2].

– Staining technique

The restoration must be free of dirt and grease. Avoid any contamination after cleaning.

a) Stain/glaze firing using IPS e.max CAD Crystall./Glaze Spray

Secure the restoration on the pin with a little IPS Object Fix Putty or Flow. Shake the IPS e.max CAD Crystall./Glaze Spray vigorously immediately before use until the mixing ball in the container is moving freely (approximately 20 seconds). If the product is not sufficiently shaken, mainly the propellant is discharged with a spraying burst. As a result, the restoration will be inadequately coated with glazing powder. Observe a distance of 10 cm between the nozzle and the surface to be sprayed. Hold the spray can as upright as possible during spraying. Spray the restoration from all sides with short bursts while simultaneously rotating the restoration so that an even, covering layer is created. Shake the can again between individual bursts of spray. Wait briefly until the glaze layer is dry and has assumed a whitish colour. Subsequently, spray the restoration a second time from all sides with short bursts while simultaneously rotating the restoration so that an even, covering layer is created. If IPS e.max Ceram Glaze Spray accidentally reaches the inner aspects of the restoration, remove it with a dry brush.

Subsequently, place the restoration on the IPS e.max CAD Crystallization Tray and fire using the stipulated parameters. Observe the directions under "How to proceed after firing".

b) Stain/glaze firing using IPS e.max CAD Crystall./Glaze Paste/Fluo

Extrude the ready-to-use glazing paste (Glaze) from the syringe and mix it thoroughly. The Glaze may be slightly thinned using IPS e.max CAD Crystall./Glaze Liquid. Apply an even layer of glazing paste with a brush, covering all areas that are to be glazed. Do not apply the Glaze too thickly. Avoid "pooling", especially on the occlusal surface of crowns. Too thin a Glaze layer may lead to an unsatisfactory gloss.

Extrude the ready-to-use shades and stains from the syringe and mix them thoroughly. The shades and stains may be slightly thinned using IPS e.max CAD Crystall./Glaze Liquid. However, the consistency should still remain pasty. Apply the mixed shades and stains directly to the unfired Glaze layer using a fine brush. More intensive shades are achieved by several staining procedures and repeated firing, not by applying thicker layers. Use IPS e.max CAD Crystall./Shades Incisal to imitate the incisal area and create the effect of translucency on the crown in the incisal and occlusal third. The cusps and fissures can be individualized using stains.

Subsequently, place the restoration on the IPS e.max CAD Crystallization Tray. Secure the restoration on the pin with a little IPS Object Fix Putty or Flow and fire using the stipulated parameters. Observe the directions under "How to proceed after firing".

c) Stain/glaze firing using IPS Ivocolor materials

Slightly wet the surface with IPS Ivocolor Mixing Liquid. This facilitates the application of the Shade pastes. Mix the Shades and Essences to the desired consistency using the respective liquids. Apply a thin coat of Shade Dentin paste to the dentin area. More intensive shades are achieved by several staining procedures and repeated firing, not by applying thicker layers. Use Shade Incisal to imitate the incisal area and create the effect of translucency on the crown in the incisal and occlusal third. The cusps and fissures can be individualized using Essences. Secure the restoration on the firing pin of the honey-comb tray with a little IPS Object Fix Putty or Flow for firing.

Mix the glazing material with the IPS Ivocolor Mixing Liquid allround or longlife to achieve the desired consistency. Thinning the glazing paste will influence its gel structure, and thus its stability and flow and application behaviour. Excessive thinning will reduce the stability of the Glaze and cause it to become runny. Moreover, if excessively thinned Glaze is applied, the degree of gloss of the fired restorations may be reduced. In an undiluted state, the paste may be difficult to apply evenly due to its gel-like consistency. However, this property of the paste can be used to selectively accentuate the restoration surface.

Apply the glazing paste in an evenly covering layer on the restoration. Secure the restoration on the firing pin of the honey-comb tray with a little IPS Object Fix Putty or Flow for firing.

How to proceed after firing

Remove the restoration from the furnace once the firing cycle is complete (wait for the acoustic signal of the furnace) and allow it to cool to room temperature in a place protected from draft. Hot objects must not be touched with metal tongs. Remove the restoration from the hardened IPS Object Fix Putty / Flow. Remove any residue with a brush and clean with water in an ultrasonic bath or with the steam jet. Residue must not be removed by blasting with Al_2O_3 or glass polishing beads. If the restoration needs to be adjusted by grinding^[2], make sure that no overheating of the ceramic occurs. Finally, polish the ground areas to a high-gloss finish.

Shade combination table

The values indicated are guiding values. When using Shades not available as blocks, please proceed as follows:

- Select the closest (brighter) block shade
- Use IPS e.max CAD Crystall./Shades and Stains or IPS Ivocolor Shades and Essences for staining or to achieve a more intensive dentin shade.

IPS e.max CAD Crystall./:

A-D shade	BL1	BL2	BL3	BL4	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4	
Shade	0				1					2					4				4		
Shade Incisal	I1, I2																				
Stains	white, creme, sunset, copper, olive, khaki, mahogany																				

IPS Ivocolor:

IPS Invocolor:

	Initial mixing										Touching up dried-out pastes																			
Components	Mixing Liquid allround					Mixing Liquid longlife					Essence Fluid					Mixing Liquid allround					Mixing Liquid longlife					Essence Fluid				
Essence	✓					✓					✓					✓					✓					-				
Shade	✓					✓					-					✓					✓					-				
Glaze Powder FLUO	✓					✓					-					✓					✓					-				
Glaze Paste FLUO	✓					✓					-					✓					✓					-				

A-D shade	BL1	BL2	BL3	BL4	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4																
Shade	SD 0				SD 1				SD 2				SD 3				SD 4				SD 5				SD 6				SD 7				SD 8			
Shade Incisal	SI 1																				SI 2				SI 3											
Essence	E 01 white, E 02 creme, E 03 lemon, E 04 sunset, E 05 copper, E 06 hazel, E 07 olive, E 08 khaki, E 09 terracotta, E 10 mahogany, E 11 cappuccino, E 12 espresso, E 13 terra, E 14 profundo, E 15 ocean, E 16 sapphire, E 17 anthracite, E 18 black, E 19 rose, E 20 coral, E 21 basic red, E 22 basic yellow, E 22 basic blue																																			

Firing parameters

Stain/glaze firing with IPS e.max CAD Crystall./Shades, Stains, Glaze

Furnaces Programat	Stand-by temperature B [°C/°F]	Closing time S [min]	Heating rate t1 [°C/°F/min]	Firing temperature T1 [°C/°F]	Holding time H1 [min]	Heating rate t2 [°C/°F/min]	Firing temperature T2 [°C/°F]	Holding time H2 [min]	Vacuum 1 11 12 [°C/°F]	Vacuum 2 21 22 [°C/°F]	Long-term cooling L [°C/°F]	Cooling rate t1 [°C/°F/min]
P300 P500 P700	403/757	6:00	90/162	820/1508	00:10	30/54	840/1544	3:00	550/820 1022/1508	820/840 1508/1544	700/1292	0
P310 G2 P510 G2 P710 G2	403/757	6:00/IRT	90/162	830/1526	00:10	30/54	850/1562	3:00	550/830 1022/1526	830/850 1526/1562	710/1310	0
CS/CS2/ CS3/CS4/ CS6	Select the corresponding program											

Staining technique with IPS Ivocolor Shades, Essence, Glaze

	Furnaces Programat	Stand-by temperature B [°C/°F]	Closing time S [min]	Heating rate t1 [°C/°F/min]	Firing temperature T1 [°C/°F]	Holding time H1 [min]	Heating rate t2 [°C/°F/min]	Firing temperature T2 [°C/°F]	Holding time H2 [min]	Vacuum 1 11 12 [°C/°F]	Vacuum 2 21 22 [°C/°F]	Long-term cooling L [°C/°F]	Cooling rate t1 [°C/°F/min]
Stain/ Glaze firing	P300 P500 P700 P310 G2 P510 G2 P710 G2	403/757	6:00	60/108	710/1310	1:00	–	–	–	450/709 842/1308	–	0	0

3 Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclar.com, and your responsible competent authority.
- The current Instructions for Use and explanation of symbols are available on the website: www.ivoclar.com/eIFU.
- The Summary of Safety and Clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>.
- Basic UDI-DI: 76152082ACERA007F4

Warnings

- IPS Natural Die Material Separator contains hexane. Hexane is highly flammable and detrimental to health. Avoid contact of the material with the skin and eyes. Do not inhale vapours and keep away from sources of ignition.
- Do not inhale ceramic dust during finishing. Use an extraction unit and wear a face mask as well as protective goggles and gloves.
- IPS e.max CAD Crystall./Glaze Spray: Do not inhale the spray mist. The container is pressurized. Protect from direct sunlight and temperatures over 50°C/122 °F. Do not forcibly open or burn the can, even after it has been completely emptied.
- If applied frequently or for a prolonged period of time, peroxide-based professional bleaching agents (carbamide peroxide; hydrogen peroxide) as well as acidulated phosphate fluorides used for caries prevention can make the surface of existing IPS e.max CAD restorations rough and matte.
- Observe the Safety Data Sheet (SDS) (available at www.ivoclar.com).

Disposal information

Remaining stocks or removed restorations must be disposed of according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks.

The following clinical residual risks are known:

- Chipping / fracture / decementation of the restorative material may lead to accidental swallowing or inhalation and dental re-treatment.
- Cement excess may lead to irritation of the soft tissue / gingiva. Progressive inflammation may lead to bone resorption and periodontitis.

4 Shelf life and storage

This product does not require any special storage conditions.

5 Additional information

Keep material out of the reach of children!

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the product for its suitability and use for any purpose not explicitly stated in the Instructions.

^[1] e.g. CEREC/inLab, PlanMill. The complete list is available at www.ivoclar.com. CEREC/inLab and PlanMill are not registered trademarks of Ivoclar Vivadent AG.

^[2] Ivoclar Flow Chart "Recommended grinding instruments for extraoral and intraoral use"

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IPS Ivocolor:

Components	Initial mixing			Touching up dried-out pastes		
	Mixing Liquid allround	Mixing Liquid longlife	Essence Fluid	Mixing Liquid allround	Mixing Liquid longlife	Essence Fluid
Essence	✓	✓	✓	✓	✓	-
Shade	✓	✓	-	✓	✓	-
Glaze Powder FLUO	✓	✓	-	✓	✓	-
Glaze Paste FLUO	✓	✓	-	✓	✓	-

A-D shade	BL1	BL2	BL3	BL4	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
Shade	SD 0				SD 1	SD 2				SD 3	SD 4	SD 5		SD 6				SD 7		SD 8
Shade Incisal	SI 1																			
Essence	E 01 white, E 02 creme, E 03 lemon, E 04 sunset, E 05 copper, E 06 hazel, E 07 olive, E 08 khaki, E 09 terracotta, E 10 mahogany, E 11 cappuccino, E 12 espresso, E 13 terra, E 14 profundo, E 15 ocean, E 16 sapphire, E 17 anthracite, E 18 black, E 19 rose, E 20 coral, E 21 basic red, E 22 basic yellow, E 22 basic blue																			

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Furnaces Programat	Stand-by temperature B [°C/F]	Closing time S [min]	Heating rate t1 [°C/min]	Firing temperature T1 [°C]	Holding time H1 [min]	Heating rate t2 [°C/min]	Firing temperature T2 [°C]	Holding time H2 [min]	Vacuum 1 T1 T2 [°C]	Vacuum 2 T1 T2 [°C]	Long-term cooling L [°C]	Cooling rate t1 [°C/min]
P300 P500 P700	403	6:00	90	820	00:10	30	840	3:00	550/820	820/840	700	0
P310 G2 P510 G2 P710 G2	403	6:00/IRT	90	830	00:10	30	850	3:00	550/830	830/850	710	0
CS/CS2/ CS3/CS4/ CS6	Select the corresponding program											

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Furnaces Programat	Stand-by temperature B [°C]	Closing time S [min]	Heating rate t1 [°C/min]	Firing temperature T1 [°C]	Holding time H1 [min]	Heating rate t2 [°C/min]	Firing temperature T2 [°C]	Holding time H2 [min]	Vacuum 1 T1 T2 [°C]	Vacuum 2 T1 T2 [°C]	Long-term cooling L [°C]	Cooling rate t1 [°C/min]
Stain/ Glaze firing	403	6:00	60	710	1:00	–	–	–	450/709	–	0	0

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