IPS Style® Ceram

[en] Instructions for Use

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1 Intended use

Intended purpose

Ceramic materials for the veneering of metal frameworks (CTE range of $13.8 - 15.2 \times 10^{-6}/K$ (25–500 °C); fabrication of layered veneers

Patient target group

- Patients with permanent teeth

Intended users

- Dentists (clinical procedure)
- Dental laboratory technicians (fabrication of restorations in the dental laboratory)

Special training

None

Use

For dental use only.

Description

IPS Style® Ceram is a feldspar-free glass-ceramic which contains leucite, fluorapatite and oxyapatite. The veneering ceramic is suitable for the fabrication of metal-ceramic restorations at firing temperatures below 880 °C. IPS Style Ceram can be used for the veneering of alloys with a CTE range of $13.8 - 15.2 \times 10^{-6}/K$ (25–500 °C) and for the fabrication of laminate veneers on refractory dies.

Product name	Product description
IPS Style Ceram Paste Opaquer IPS Style Ceram Intensive Paste Opaquer IPS Style Ceram Powder Opaquer 870 IPS Style Ceram Intensive Powder Opaquer 870	Masking of metal frameworks and adjustment of the base shade
IPS Style Ceram Margin	Fabrication of ceramic margins on metal frameworks
IPS Style Ceram Deep Dentin IPS Style Ceram Dentin IPS Style Ceram Incisal	Veneering of metal frameworks; fabrication of layered veneers in the anterior region
IPS Style Ceram Transpa IPS Style Ceram Mamelon IPS Style Ceram Occlusal Dentin IPS Style Ceram Incisal Edge IPS Style Ceram Opal Effect IPS Style Ceram Special Incisal IPS Style Ceram Inter Incisal IPS Style Ceram Cervical Transpa	Characterization of ceramic veneers on metal frameworks; characterization of layered veneers in anterior teeth
IPS Style Ceram One	Veneering of metal frameworks
IPS Style Ceram Add-On Margin IPS Style Ceram Add-On Dentin IPS Style Ceram Add-On Incisal IPS Style Ceram Add-On Bleach IPS Style Ceram Add-On 690°C	Adjustment of ceramic margins and restorations on metal frameworks
IPS Style Ceram Gingiva IPS Style Ceram Intensive Gingiva IPS Style Ceram Basic Gingiva	Veneering of metal frameworks in the gingival area Characterization of Gingiva materials

Accessories for the medical device

Product name	Product description
IPS Paste Opaquer Liquid	For diluting IPS Style Ceram Paste Opaquer
IPS Powder Opaquer Liquid	For mixing IPS Style Ceram Powder Opaquer
IPS Margin Build-Up Liquid allround IPS Margin Build-Up Liquid carving	For mixing IPS Style Ceram Margin layering materials
IPS Build-Up Liquid (allround)	For mixing IPS Style Ceram layering materials

Indications

- Missing tooth structure in anterior and posterior teeth
- Partial edentulism in the anterior and posterior region

Types of restorations:

- Crowns
- 3–14-unit bridges
- Veneers

Contraindications

- The use of this product is contraindicated if the patient is known to be allergic to any of its ingredients.

Limitations of use

- Untreated bruxism (a splint is indicated after incorporation)
- Mixing and processing together with other dental ceramics
- Failure to observe the necessary minimum connector and framework thicknesses
- Exceeding or falling short of the stipulated veneering layer thicknesses
- Failure to observe the layer thickness ratio between the framework and layering ceramic
- Veneering of dental alloys not within the stipulated CTE range
- Veneering of titanium and zirconium oxide frameworks
- The final restoration must not be reused.

Side effects

There are no known side effects to date.

Interactions

There are no known interactions to date.

Clinical benefit

- Reconstruction of chewing function
- Restoration of esthetics

Composition

	Product components powders and pastes	Composition
Opaquer	IPS Style Ceram Paste Opaquer IPS Style Ceram Intensive Paste Opaquer	- Mixed glass-ceramics containing leucite and oxyapatite crystals - Glass - Zirconium oxide - Inorganic pigments - Glycols (butanediol, glycerine) - Polyvinylpyrrolidone - Highly dispersed silicic acid
	IPS Style Ceram Powder Opaquer 870 IPS Style Ceram Intensive Powder Opaquer 870	- Mixed glass-ceramics containing leucite and oxyapatite crystals - Zirconium oxide - Inorganic pigments
Margin	IPS Style Ceram Margin	- Mixed glass-ceramics containing leucite, fluorapatite and oxyapatite crystals - Inorganic pigments
Dentin	IPS Style Ceram Deep Dentin	- Mixed glass-ceramics containing leucite and oxyapatite crystals
	IPS Style Ceram Dentin	- Mixed glass-ceramics containing leucite, fluorapatite and oxyapatite crystals - Inorganic pigments
Incisal	IPS Style Ceram Incisal	- Mixed glass-ceramics containing leucite and oxyapatite crystals - Inorganic pigments
Impulse	IPS Style Ceram Occlusal Dentin IPS Style Ceram Mamelon IPS Style Ceram Inter Incisal IPS Style Ceram Incisal Edge IPS Style Ceram Gingiva IPS Style Ceram Basic Gingiva IPS Style Ceram Intensive Gingiva	- Mixed glass-ceramics containing leucite and oxyapatite crystals - Mixed glass-ceramics containing leucite, fluorapatite and oxyapatite crystals - Inorganic pigments
	IPS Style Ceram Transpa IPS Style Ceram Special Incisal IPS Style Ceram Incisal Edge IPS Style Ceram Cervical Transpa	- Mixed glass-ceramics containing leucite and oxyapatite crystals - Inorganic pigments

	Product components powders and pastes	Composition
Effect	IPS Style Ceram Opal Effect	– Mixed glass-ceramics containing leucite and oxyapatite crystals – Mixed glass-ceramics containing leucite, fluorapatite and oxyapatite crystals – Leucite and fluorapatite glass-ceramics – Inorganic pigments
One	IPS Style Ceram One	– Mixed glass-ceramics containing leucite and oxyapatite crystals – Mixed glass-, -ceramics containing leucite, fluorapatite and oxyapatite crystals – Inorganic pigments
Add-on	IPS Style Ceram Add-On Margin	– Glass – Leucite glass-ceramic – Mixed glass-ceramics containing leucite, fluorapatite and oxyapatite crystals – Inorganic pigments
	IPS Style Ceram Add-On Dentin	
	IPS Style Ceram Add-On Bleach	
	IPS Style Ceram Add-On Incisal	
	IPS Style Ceram Add-On 690 °C	– Glass – Leucite glass-ceramic – Inorganic pigments

Accessories for the medical device

Product components ceramic mixing liquids	Composition
IPS Paste Opaquer Liquid	– Glycol (butanediol)
IPS Powder Opaquer Liquid	– Water – Glycol (propanediol) – Carbonic acid – Metallic salt
IPS Build-Up Liquid (allround)	– Water – Glycol (butanediol) – Metallic salt
IPS Margin Build-Up Liquid (allround)	– Water – Cellulose derivative – Metallic salt
IPS Margin Build-Up Liquid (carving)	

2 Application

! NOTICE

Observe the following specifications:

- Wall thickness for metal frameworks: Crowns min. 0.3 mm; abutment crowns min. 0.5 mm.
- Ceramic layering thickness: min. 0.8 mm; max. 1.7 mm.

! NOTICE

Only use the system components mentioned in chapter Intended use. Proceed as follows:

- Mix powder materials with mixing liquid before use.
- Extrude and mix paste materials before use. Use mixing liquid for diluting.

! NOTICE

Conduct all firing cycles according to the instructions of the manufacturer. Firing parameters can be found in chapter Technical specifications.

2.1 Layering crowns/bridges

Fabricating the framework

1. Fabricate and oxidize the framework with a reduced anatomical shape according to the instructions of the alloy manufacturer.

Applying and firing the opaquer

1. Apply a thin, non-covering layer of paste or powder opaquer on the framework.
2. Conduct first Opaquer firing (Wash firing).
3. Apply a thin, covering layer of paste or powder opaquer on the framework.
4. Apply an even covering layer of Gingiva paste or powder opaquer on the gingiva portions.

5. Customize the framework using Intensive paste or powder opaquer.
6. Conduct second Opaquer firing.

Layering and firing Gingiva material

1. Layer the gingiva area.
2. Optional: Customize the restoration using characterization material.
3. Conduct Gingiva firing.
4. Repeat the steps 1 to 3.

A One-layer technique

Layering and firing One material

1. Fully build up the restoration using the one-layer material
2. Conduct One firing.
3. Repeat the steps 1 to 2.

B Standard and individual layering technique

In case of ceramic margins: Layering and firing the ceramic margin

1. Isolate and seal the die in the margin areas.
2. Apply Margin material in the cervical area.
3. Dry the Margin material.
4. Lift the framework with Margin material off the die.
5. Conduct Margin firing.
6. Check the fit.
7. Remove interfering areas.
8. Repeat the steps 2 to 7.

Layering and firing Dentin and Incisal material

1. If space is limited: Apply opaque Dentin material on the opaquerized framework.
2. Completely layer the restoration with Dentin material.
3. Cut back the restoration.
4. Add Incisal material to the incisal area.
5. Optional: Customize the restoration using characterization material.
6. Conduct Dentin and Incisal firing.
7. Repeat the steps 4 to 6.

Characterizing, glazing and firing the restoration with IPS Ivocolor

NOTICE! Make sure that the adhesive surface of the restoration is completely free of ceramic material.

1. Finalize the anatomical shape and surface structure.
2. Customize the restoration using stains.
3. Conduct Stains firing.
4. Apply glaze material using a brush.
5. Conduct Glaze firing.

Applying and firing corrections

NOTICE! Corrections can be fired with the glaze firing cycle or after the glaze firing cycle.

1. **Correction with Glaze firing:** Apply the Add-On material.
2. Conduct Glaze firing.
3. **Correction after glaze firing:** Apply the Add-On 690° material.
4. Conduct corrective firing.

2.2 Layering veneers on refractory dies

Duplicating the die

1. Duplicate the die in silicone.
2. Remove the die.
3. Pour a commercial refractory die material into the silicone mould.
4. Degas the refractory die.
5. Soak the die in water.

Applying and firing the Wash

1. Mix IPS Style Ceram Deep Dentin, Dentin, Incisal or Transpa layering material with mixing liquid.
2. Apply a thin layer of the material on the refractory die.
3. Conduct Wash firing.

Layering and firing the veneer

1. Layer Dentin and Impulse material on the die.
2. Conduct Dentin/Impulse firing.
3. Build up incisal areas using Incisal materials.
4. Conduct Incisal firing.

Characterizing, glazing and firing the restoration

See chapter Layering crowns/bridges.

Divesting the veneer

1. Remove the die material from the veneer using a polishing jet medium (max. 1 bar).

Conditioning the veneer for adhesive cementation

1. Blast the inner aspect of the veneer with Al_2O_3 (max. 1 bar).
2. Clean the veneer with the steam jet.
3. Dry the veneer with oil-free air.
4. Etch the inner aspect of the veneer with IPS Ceramic Etching Gel for 60 s.
5. Repeat the steps 2 to 3.

Seating the restoration

NOTICE! Seating is exclusively done by the dentist.

- Seat crowns and bridges using conventional or adhesive methods.
- Seat veneers using adhesive methods.

3 Technical specifications

Firing table

NOTICE! These firing parameters are guidance values. They are valid for the Programat® furnaces from Ivoclar Vivadent AG.

Deviations (approx. +/- 10 °C) may occur:

- depending on the furnace generation
- when using ceramic furnaces from other manufacturers
- due to regional differences in the power supply
- if several electrical devices are operated on the same circuit

IPS Style® Ceram – Firing parameters

IPS Style® Ceram	Firing temperature T [°C]	Stand-by temperature B [°C]	Closing time S [min]	Heating rate t [°C/min]	Holding time H [min]	Vacuum on V1 [°C]	Vacuum off V2 [°C]
1 st /2 nd Powder opaquer firing	870	403	4:00	100	1:00	450	869
1 st /2 nd Paste opaquer firing	870	403	6:00	100	2:00	450	869
1 st /2 nd Margin firing	840	403	6:00	60	1:30	450	839
1 st Dentin / Incisal / Gingiva / One firing	790	403	6:00	60	1:00	450	789
2 nd Dentin / Incisal / Gingiva / One firing	780	403	6:00	60	1:00	450	779
Stains firing using IPS Ivocolor*	750	403	6:00	60	1:00	450	749
Glaze firing / Add-On with Glaze firing	750	403	6:00	60	1:00	450	749
Add-On after Glaze firing	690	403	6:00	60	1:00	450	689

* Stains firing is used to fixate IPS Ivocolor Shades/Essences in place and is particularly recommended for comprehensive characterizations (e.g. one-layer technique).

IPS Style® Ceram Veneer – Firing parameters

IPS Style® Ceram Veneer	Firing temperature T [°C]	Stand-by temperature B [°C]	Closing time S [min]	Heating rate t [°C/min]	Holding time H [min]	Vacuum on V1 [°C]	Vacuum off V2 [°C]
Wash firing	810	403	8:00	50	1:00	450	809
Cervical firing	800	403	8:00	50	1:00	450	799
Dentin / Impulse firing	800	403	8:00	50	1:00	450	799
Incisal firing	800	403	8:00	50	1:00	450	799
Stain/glaze firing using IPS Ivocolor materials	750	403	8:00	50	1:00 – 1:30	450	749

Properties of IPS Style Ceram according to ISO 6872:2015 + Amd 1:2018 and ISO 9693:2019. IPS Style Ceram has been classified as a Type I, Class 1 dental ceramic.

Physico-chemical properties of the product components of IPS Style Ceram according to 1) ISO 6872:2015 + Amd 1:2018 and 2) ISO 9693:2019.

Product components/ properties	Radioactivity [Bq/g U ²³⁸]	Biaxial flexural strength σ [MPa]	Chem. solubility [$\mu\text{g}/\text{cm}^2$]	Adhesive bond Tb [MPa]
	ISO 6872:2015 + Amd 1:2018			ISO 9693:2019
IPS Style Ceram Paste Opaquer	≤ 1.0	≥ 50	< 100	> 25
IPS Style Ceram Intensive Paste Opaquer				
IPS Style Ceram Powder Opaquer 870				
IPS Style Ceram Intensive Powder Opaquer 870				
IPS Style Ceram Deep Dentin				
IPS Style Ceram Dentin				
IPS Style Ceram Occlusal Dentin				
IPS Style Ceram Mamelon				
IPS Style Ceram Incisal				
IPS Style Ceram Incisal Edge				
IPS Style Ceram Inter Incisal				
IPS Style Ceram Special Incisal				
IPS Style Ceram Transpa				
IPS Style Ceram Cervical Transpa				
IPS Style Ceram Gingiva				
IPS Style Ceram Intensive Gingiva				
IPS Style Ceram Basic Gingiva				
IPS Style Ceram One				
IPS Style Ceram Margin				n.a.
IPS Style Ceram Opal Effect				
IPS Style Ceram Add-On Margin				
IPS Style Ceram Add-On Dentin				
IPS Style Ceram Add-On Incisal				
IPS Style Ceram Add-On BL				
IPS Style Ceram Add-On 690°C				

Thermal properties of the product components of IPS Style Ceram according to 1) ISO 6872:2015 + Amd 1:2018 and 2) ISO 9693:2019.

Product components/properties	Coefficient of linear thermal expansion (CTE)				Glass transition temperature t_g [°C]	Thermocycling test
	2 firing cycles [10 ⁻⁶ /K]	4 firing cycles [10 ⁻⁶ /K]	Average value [10 ⁻⁶ /K]	Measuring range [°C]		
						ISO 9693:2019
IPS Style Ceram Paste Opaquer	13.5 ± 0.5	13.7 ± 0.5	13.6 ± 0.5	25 - 500	540 ± 20	applicable
IPS Style Ceram Intensive Paste Opaquer						
IPS Style Ceram Powder Opaquer 870						
IPS Style Ceram Intensive Powder Opaquer 870						
IPS Style Ceram Deep Dentin	12.0 ± 0.5	12.4 ± 0.5	12.2 ± 0.5		525 ± 20	applicable ¹ / n.a. ¹
IPS Style Ceram Dentin						
IPS Style Ceram Occlusal Dentin						
IPS Style Ceram Mamelon						
IPS Style Ceram Incisal						
IPS Style Ceram Incisal Edge						
IPS Style Ceram Inter Incisal						
IPS Style Ceram Special Incisal						
IPS Style Ceram Transpa						
IPS Style Ceram Cervical Transpa						
IPS Style Ceram Gingiva						
IPS Style Ceram Intensive Gingiva						
IPS Style Ceram Basic Gingiva						
IPS Style Ceram One						
IPS Style Ceram Margin	12.8 ± 0.5	13.1 ± 0.5	13.0 ± 0.5	580 ± 20	n.a.	
IPS Style Ceram Opal Effect	12.0 ± 0.5	12.4 ± 0.5	12.2 ± 0.5	510 ± 20		
IPS Style Ceram Add-On Margin	12.7 ± 0.5	13.3 ± 0.5	13.0 ± 0.5	25 - 400		480 ± 20
IPS Style Ceram Add-On Dentin	13.0 ± 0.5	13.6 ± 0.5	13.3 ± 0.5	475 ± 20		
IPS Style Ceram Add-On Incisal						
IPS Style Ceram Add-On BL						
IPS Style Ceram Add-On 690°C	12.0 ± 0.5	12.4 ± 0.5	12.2 ± 0.5	455 ± 20		

¹ Applicable to IPS Style Ceram Dentin, IPS Style Ceram Incisal and IPS Style Ceram Transpa; all others: n.a.

4 Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, www.ivoclar.com and your responsible competent authority.
- The current Instructions for Use and the explanation of symbols are available on the website: www.ivoclar.com/eIFU
- The Summary of Safety and Clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>.
- Basic UDI-DI: 76152082ACERA004EW

Warnings

- Observe the Safety Data Sheet (SDS) (available at www.ivoclar.com).
- Do not inhale ceramic dust during finishing. Use suction equipment, a face mask and protective goggles.
- If applied frequently or for a prolonged period of time, peroxide-based professional bleaching agents (carbamide peroxide; hydrogen peroxide) as well as acidulated phosphate fluorides used for caries prevention can make the surface of existing IPS Style Ceram restorations rough and matte.
- The product contains hexane. Hexane is highly flammable and detrimental to health. Avoid contact of the material with the skin and eyes. Do not inhale vapours. Keep material away from sources of ignition.

Disposal information

Remaining stocks or removed restorations must be disposed of according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks.

The following known residual clinical risks exist:

- Chipping / fracture / debonding of the restorative material may lead to accidental swallowing or inhalation and dental retreatment

5 Shelf life and storage

- Storage temperature for ceramic pastes and mixing liquids: 2–28 °C (36–82 °F).
- Storage temperature for ceramic powders: No special storage conditions are required.
- Expiry date: See note on bottle, jar or packaging.
- Do not use the product after the indicated expiry date.
- Store in a vibration-free environment.
- Keep dry.
- Keep out of sunlight.
- Check to make sure that the packaging and the product are intact and undamaged before use. If in doubt, contact Ivoclar Vivadent AG or your local sales partner.

6 Additional information

Keep material out of the reach of children!

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

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